

Jackson-Pratt® Drains Perfecting the Art of Performance

Hemaduct® Cardiothoracic Drain and Multiport Connector

The Jackson-Pratt® brand Hemaduct® cardiothoracic drain and multiport connector offer greater patient comfort and promote earlier patient mobility over standard thoracotomy tubes.

Use the Hemaduct cardiothoracic drain and multiport connector for procedures such as CABG or OPCABG, valve replacement or repair, thoracotomy and lobectomy.

The Hemaduct cardiothoracic drain features:

- 19 Fr along the entire length of drain and tube
- Effective distribution of suction and flow over entire drain
- Innovative design providing pathways around obstructions
- Silicone material for strength and softness

The multiport connector offers:

- Connection for up to three wound drains (with caps covering unused ports)
- Connection of the Hemaduct CT drain(s) directly to standard chest drain unit reservoirs



Hemaduct cardiothoracic drain and multiport connector.



The multiport connector allows connection for up to three wound drains.



Unique system of channels, lumens and portals provides suction throughout entire drain.


Procedural Flexibility and Superior Suction

Compare the Effectiveness:

Hemaduct® Drain vs. Johnson & Johnson Blake® Drain

Independent testing shows that over a 10-min. period, Hemaduct's multiple pathways and superior suction throughout the entire length of the drain allow it to suction 10 times more fluid than the competitive channel drain.

B=Blake drain H=Hemaduct drain



Cat. No.	Description
CT-HU1900	Jackson-Pratt® brand Hemaduct® cardiothoracic drain, 19 Fr
CT-HU1910	Jackson-Pratt brand Hemaduct cardiothoracic drain with trocar, 19 Fr
CT-131	Jackson-Pratt brand Hemaduct cardiothoracic multiport connector

For more information or to place an order, please contact your Cardinal Health, Medical Products and Services sales representative or call customer service at 1.800.964.5227.



Multiport connector attaches directly to any chest drain unit.



Hemaduct cardiothoracic drain attaches directly to Jackson-Pratt 400cc bulb reservoir.

Hemaduct, Jackson-Pratt and Medi-Vac are registered trademarks of Cardinal Health, Inc. or one of its subsidiaries. Blake is a registered trademark of Johnson & Johnson.

© Copyright 2003. Cardinal Health, Inc. or one of its subsidiaries. All rights reserved.
Lit. No. MV00541 (0803/5M/MC1565)

♻️ Printed on recycled paper with 10% post-consumer fiber.

Cardinal Health
Medi-Vac® Products
1500 Waukegan Road
McGaw Park, IL 60085

www.cardinal.com

